



Auto Chart Projector **LED**
ACP·900

CE



ELLEGI MEDICAL OPTICS

ACP·900



SPECIFICHE

Test	40
Velocità selezione	0.2 sec.
Maschere	1 Aperta 8 Linee verticali 5 Linee orizzontali 21 Lettere singole
Filtri	Rosso / Verde
Distanza proiezione	da 2.0 a 7.0mt (standard: 3mt)
Ingrandimenti	30X (a 5mt)
Inclinazione	+ / - 15°
Lampadina	LED
Spegnimento automatico	Da 5 a 10 minuti
Consumo	45VA
Alimentazione	AC 100 - 240v, 50 / 60Hz
Dimensioni/peso corpo	290x186x198 mm, 3.6 Kg
Dimensioni peso telecomando	64x195x20, 110g
Interfaccia	Connettore 8 Pin per la connessione simultanea di ottotipi,
maschera e filtro R/V, tramite unità di controllo elettrica del forottero	
Accessori	1 Foderina 1 Schermo metallico polarizzato 2 Fusibili e 2 Batterie
Accessori opzionali	Lente polarizzata Occhiali rosso / verdi

Proiezione chiara e precisa

Quattro differenti tipi di grafici, immagini chiare e nitide per un test di precisione

Nuovo design

Compatto ed innovativo

TIPI DI GRAFICI (A/B/C/M)

A TYPE	
LETTERS	Letters 1: K 0.05 Letters 2: HRN 0.10, FDA 0.15 Letters 3: ZTCB 0.20, KENVH 0.30, OLFA 0.40 Letters 4: KCAEH 0.50, ZTCV 0.60, OLFA 0.70 Letters 5: TOEPL 0.80, VZTA 0.90, HAFND 1.00 Letters 6: KENV 1.20, VZTA 1.30, ZKVP 2.00
SNELLEN	Snelten 1: MEW 0.10, ME 0.15 Snelten 2: MEW 0.20, ME 0.30, ME 0.40 Snelten 3: MEW 0.50, ME 0.60, ME 0.70 Snelten 4: MEW 0.80, ME 0.90, ME 1.00
CHILDREN	Children 1: Butterfly 0.10 Children 2: Butterfly 0.20, Butterfly 0.40 Children 3: Butterfly 0.60, Butterfly 0.80, Butterfly 1.00
NUMBER	Number 1: 235 0.10, 8 0.15 Number 2: 235 0.20, 8 0.30, 8 0.40 Number 3: 235 0.50, 8 0.60, 8 0.70 Number 4: 235 0.80, 8 0.90, 8 1.00 Number 5: 235 1.20, 8 1.30, 8 2.00
FUNCTION	Red / Green Duochrome Balance Worth for Dot Schobor Cross Grid Cross Cylinder Dots Alignment Clock Dial Fixation Binocular Balance Phoria Phoria with Fixation Coincidence Vertical Coincidence Horizontal Minute Stereo

B TYPE	
LETTERS	Letters 1: E 0.08 Letters 2: HRN 0.10, FDA 0.15 Letters 3: ZTCB 0.20, KENVH 0.30, OLFA 0.40 Letters 4: KCAEH 0.50, ZTCV 0.60, OLFA 0.70 Letters 5: TOEPL 0.80, VZTA 0.90, HAFND 1.00
SNELLEN	Snelten 1: MEW 0.10, ME 0.15 Snelten 2: MEW 0.20, ME 0.30, ME 0.40 Snelten 3: MEW 0.50, ME 0.60, ME 0.70 Snelten 4: MEW 0.80, ME 0.90, ME 1.00
CHILDREN	Children 1: Butterfly 0.10 Children 2: Butterfly 0.20, Butterfly 0.40 Children 3: Butterfly 0.60, Butterfly 0.80, Butterfly 1.00
LANDOLT	Landolt 1: 0.20, 0.30, 0.40, 0.50, 0.60, 0.70, 0.80, 0.90, 1.00 Landolt 2: 0.20, 0.30, 0.40, 0.50, 0.60, 0.70, 0.80, 0.90, 1.00 Landolt 3: 0.20, 0.30, 0.40, 0.50, 0.60, 0.70, 0.80, 0.90, 1.00 Landolt 4: 0.20, 0.30, 0.40, 0.50, 0.60, 0.70, 0.80, 0.90, 1.00
FUNCTION	Red / Green Duochrome Balance Worth for Dot Schobor Cross Grid Cross Cylinder Dots Alignment Clock Dial Fixation Binocular Balance Phoria Phoria with Fixation Coincidence Vertical Coincidence Horizontal Minute Stereo

C TYPE	
LETTERS	Letters 1: E 400 Letters 2: HRN 200, FDA 150 Letters 3: AKN 100, ECD 80 Letters 4: ZTCB 70, KCAEH 40, ZTCV 30, OLFA 20 Letters 5: TOEPL 20, VZTA 15, HAFND 10
SNELLEN	Snelten 1: MEW 100, ME 50 Snelten 2: MEW 60, ME 30, ME 20 Snelten 3: MEW 40, ME 20, ME 15 Snelten 4: MEW 20, ME 15, ME 10
CHILDREN	Children 1: Butterfly 100 Children 2: Butterfly 60, Butterfly 30, Butterfly 20 Children 3: Butterfly 40, Butterfly 20, Butterfly 15
NUMBER	Number 1: 235 200, 8 150 Number 2: 235 100, 8 80, 8 30 Number 3: 235 70, 8 40, 8 20 Number 4: 235 40, 8 20, 8 15 Number 5: 235 20, 8 15, 8 10
FUNCTION	Red / Green Duochrome Balance Worth for Dot Schobor Fixation Cross Grid Cross Cylinder Dots Alignment Clock Dial Phoria Phoria with Fixation Coincidence Vertical Minute Stereo Vertical Line Horizontal Line

M TYPE	
	USA / CANADA

ELLEGI MEDICAL OPTICS S.R.L.

Via Pisciarelli, 79 - c/o Complesso Altamira - 80078 Pozzuoli (NA)
tel.: +39 081 549 96 45 - fax: +39 081 549 80 44 - info@ellegimedical.it - www.ellegimedical.it

P.IVA e Reg. Imp. Napoli 05354730631 - REA NA 438621 - Cod. Mecc. NA 044696
Capitale Sociale: 26.000,00 euro int. vers. - IBAN: IT84 D030 3203 4000 1000 0012 529